

iPulse™ IAB

Abiomed's iPulse™ intra-aortic Balloon Catheter is a 40cc, 8F, intra-aortic counterpulsation device that passively unloads the left ventricle, reduces the workload of the heart and enhances cardiac output by approximately 0.8 L/min.¹ The iPulse™ Catheter can be used on the iPulse™ Console or on other OEM consoles.



SIMPLE

- Proprietary balloon material, AngioFlex® II has a smooth surface making it simple to insert.
- Sheath or sheathless insertion.

SENSIBLE


- Enhance ventricular unloading.
- Increase cardiac output with up to 0.8 L of flow per minute.
- Compatible with other OEM consoles.
- Fastest inflate and deflate (data on file).
- More accurate trigger (data on file).

SMART

- Automatic timing.
- Automatic arrhythmia detection.
- Automatic purging and filling of balloons.



Recovering hearts. Saving lives.™



TZJ c'XVäHéZXMi 'dch

Catheter size	8F.
Catheter Length	720 mm. Insertable length
Balloon Size	40 cc.
Guidewire Size	.025 in.

Ordering Information

Part Number	Description
0056-0040	iPulse IAB

Instructions for Use (IFU) can be found on www.abiomed.com



Abiomed, Inc.
 22 Cherry Hill Drive
 Danvers, Massachusetts 01923 USA
 Voice: 978-777-5410
 Facsimile: 978-777-8411
 Email: marketing@abiomed.com

Abiomed Europe GmbH
 Neuenhofer Weg 3
 52074 Aachen, Germany
 Voice: +49 (241) 8860-0
 Facsimile: +49 (241) 8860-111
 Email: europe@abiomed.com